Orthodontic Archwires

Intended Use

Wires provide force to move the teeth into the desired Class I bite and to maintain the formation of the teeth until the roots of the teeth remain in the preferred position.

Description of the device



An orthodontic archwire is a metal wire available in different sizes and dental arch shapes and materials. Some archwires have a coating to make them aesthetic. Archwires are a part of an orthodontic appliance consisting of bands and brackets.

1. Part Numbers

Product	Part	Variants	Material	Composition
Stainless Steel Soft Wire Spool	SSSOF*	Wire Size: .016"; Length: 15'	Stainless Steel 316	Chromium 16-18%, Nickel 10-14%, Molybdenum 2-3%, Manganese <2%, Iron Balance
Stranded Stainless Steel Spool	SST*	Size and strand quantity	Stainless Steel	Aluminum <3.5%, Chromium <30%, Cobalt <30%, Copper <34%, Iron <85%, Manganese <10%, Molybdenum <18%, Nickel <35%, Silicon <4.5%, Tantalum <5.5%, Tungsten <6.5%, Vanadium <0.5%
Stranded Stainless Steel Straight Length	MS*	Size and strand quantity		
Australian Stainless Steel Archwires	AW*	Size and archform	Stainless Steel 302	Chromium 16-18%, Nickel 8-10%, Manganese <2%, Cobalt <1%, Iron Balance
Stranded Stainless Steel Archwires	C*	Size, shape and archform		
Stainless Steel Archwires	SS* S3* XR*	Size, shape and archform	Stainless Steel 304	Iron 69.5%; Chromium 18.5%; Nickel 9.0%; Manganese 1.0%; Cobalt 0.75%; <1% Silicon, Carbon, Phosphorus, Sulfur
Posted Stainless Steel Archwires	PS* PX*	Size, archform and post placement	Stainless Steel 304	Iron 69.5%; Chromium 18.5%; Nickel 9.0%; Manganese 1.0%; Cobalt 0.75%; <1% Silicon, Carbon, Phosphorus, Sulfur
Looped Stainless Steel Archwires Looped Stainless Steel Archwires	KH* TL*	Size, archform loop type, and loop position		
Stainless Steel Straight Lengths	STSS*	Size and shape		
TitanMoly/Beta Titanium Archwires	TM* B3*	Size, shape, archform, and finish	Titanium	Molybdenum 10-20%, Zirconium 4.5-10%, Tin 3.75-8%, Titanium Balance
TitanMoly/Beta Titanium Looped Archwires	TMTL* B3TL*	Size, archform: and Loop Placement		
TitanMoly/Beta Titanium Straight Lengths	STB* STTM*	Size and shape		
Colboloy Blue Archwires	CB*	Size, shape, and archform	Cobalt	Cobalt 39-41% Chromium 19-21%

Colboly Blue Straight Lengths	STRCB*	Size and shape		Nickel 15-16% Iron 12-16% Molybdenum 6.5-7.5% Manganese 1.5-2% Silicon <1.2% Carbon <0.2
Nickel Titanium Archwires	G4* SE*	Size, shape, archform and finish	Nickel Titanium	Nickel 54.5 – 57%; Carbon < 0.05%; Cobalt <0.05%; Copper <0.05%; Chromium <0.05%; Hydrogen 0.05%; Iron 0.05%; Niobium <0.025%; Nitrogen plus Oxygen <0.05%; Any Single Trace Element <0.1%; Total Trace Elements <0.25%; Titanium balance%
Nickel Titanium Thermal Archwires	T* CU* M5*	Size, shape, archform and finish		
Nickel Titanium Pretorque	PT*	Size, archform: torque placement and finish:		
Nickel Titanium Reverse Curve	RC*	Shape and archform		
Nickel Titanium Straight Lengths	SENT*	Size and shape		
Nickel Titanium Spooled	SESP*	Size		

Warnings and Precautions

General		Devices are designed and manufactured for single patient use only. The reuse or reprocessing may cause cross infection or loss of mechanical properties due to natural wear of the product.
General		Appliance breakage and or distortion require immediate attention by the orthodontist.
General		Patients should wear a mouthguard when engaging in sports activities; injury to the patient and /or breakage or distortion of the appliance can occur. In the event of injury or ingestion, instruct the patient to seek medical care.
General		Care should be taken to avoid habits or foods that could break or distort the wire. Follow orthodontist recommendations for proper at home care. In the event of injury or ingestion, instruct the patient to seek medical care.
General		This product contains nickel or chromium, a chemical known to the State of California to cause cancer. In the event of an allergic reaction, immediately remove orthodontic appliance and instruct the patient to seek medical care.
Coated Archwires		The coating on the archwire may chip during regular use.
Stainless Steel		This medical device is MR Conditional. Scanning of the head and neck may cause image artifacts. If an MRI is needed, contact the radiologist prior to the procedure for instructions regarding appliance removal, if required.

Clinical Benefits

Archwires, in conjunction with orthodontic appliances, brings an aesthetic improvement with the alignment and leveling of the teeth, better hygiene, more efficient masticatory function facilitating swallowing, better digestion, prevention of occlusal problems, better speech and breathing.

Residual Risks

Orthodontic treatments are subject to relapse. Archwires, as part of the fixed orthodontic system, can

contribute to discoloration or decalcification of teeth, bone and root resorption, and periodontal complications

Potential Side Effects or Adverse Events

- Allergic reaction, particularly if wire contains nickel.
- Discomfort and pain
- Oral and mucosal damage
- Ingestion of a section of broken wire
- Adverse tooth movement resulting from wire distortion
- Archwires, as part of the fixed orthodontic system, can contribute to difficulties in maintaining oral hygiene, difficulty speaking or chewing

Instructions for use

The use of Orthodontic wire is restricted to dental professionals with qualification in orthodontics. Use by persons outside this context may result in irreversible damage to the patient. The choice of material to be used is made by the qualified professional, based on their knowledge acquired during their training in conjunction with the diagnosis and treatment plan for each patient.

1. Select the Wire to be used for the particular arch.
2. Place the archwire into the bracket slot. Reference archwire midline markings where applicable.
3. Tie the wire with elastomeric or ligature tie (not applicable for self-ligating brackets)
4. Cut off the excess wire at the distal of the furthest most attachment, or turn the end of the wire 90 degrees to the lingual (heat treatment optional according to orthodontist training), to prevent injuries in soft tissue of the patient.
5. Patient should be monitored regularly when undergoing orthodontic treatment. Particular care should be taken to monitor for wire distortion.

Any serious incident experienced with the use of the product should be reported to the manufacture and the competent authority.

Maintenance and Storage instructions

Store archwire in product packaging until use to avoid wire distortion.

Disposal

The professional should use good clinical practices for the disposal of medical products in accordance with current sanitary regulations.

Patient Instructions

	Patients should wear a mouthguard when engaging in sports activities; injury to the patient and /or breakage or distortion of the appliance can occur. In the event of injury or ingestion, instruct the patient to seek medical care.
	Practice good oral hygiene to prevent plaque installation and subsequent periodontal disease.
	Appliance breakage and or distortion requires immediate attention by the orthodontist.
	Care should be taken to avoid habits or foods that could break or distort the wire. Follow orthodontist recommendations for proper at home care. In the event of injury or ingestion, instruct the patient to seek medical care.
	This product contains nickel or chromium, a chemical known to the State of California to cause cancer. In the event of an allergic reaction, immediately remove orthodontic appliance and instruct the patient to seek medical care.

This medical device is MR Conditional. Scanning of the head and neck may cause image artifacts. If an MRI is needed, contact the radiologist prior to the procedure for instructions regarding appliance removal, if required.

Warranty

Limited Warranty. Product(s) manufactured by G&H are warranted to be free of defects in materials and manufacture for ninety (90) days from date of delivery. Any material or manufacture defect covered by this limited warranty which occurs during normal use and is reported to Seller in writing during the period of ninety (90) days from the date such products are shipped to customer. G&H's obligation hereunder, upon verification of the defect or error, shall be to provide one of the following, in G&H's sole discretion: (i) replacement at no charge to customer; or (ii) credit the purchase price to customer. If, upon the inspection of G&H product to which this warranty applies, G&H determines that a claimed defect was not due to its manufacture or materials, G&H will invoice customer for the value of product and services at customer's expense. This warranty shall not apply to any product that has been subject to misuse, misapplication, neglect (including but not limited to use of unauthorized parts or attachments), adjustment or repair and/or improper handling and storage. Identifiable items manufactured by others but installed in or affixed to G&H products, or distributed by G&H to customer, are not warranted by G&H but bear only the express warranty, if any, of the manufacturers thereof. No representative of G&H has any authority to waive, alter, vary or add to the terms hereof without prior approval in writing, to customer, signed by an officer of G&H. G&H's liability for its products, whether for breach of contract, negligence, strict liability in tort, or otherwise, shall be limited, in G&H's sole discretion, to the replacement of the products or parts thereof or a credit to the purchase price to customer. G&H will not be liable for any other injury, loss, damage or expense, whether direct or consequential, including but not limited to loss of use, income, profit or production, or increased cost of operation, or spoilage of or damage to material, arising in connection with the sale, installation, use of, inability to use, or the repair or replacement of, or late delivery of, G&H products.

Limited Warranty. Product(s) manufactured by Highland Metals are warranted to be free from defects. Products must be returned within 90 days from the date of purchase and must be returned in original packaging and in resalable condition. HMI™ reserves the right to final determination as to acceptance of product exchange or credit. Restock fee may apply. Prior to the return of any merchandise, please contact HMI™ Customer Service for a return authorization number

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