**Orthodontic Bands**

**Intended Use**

Orthodontic bands provide a base to affix weldable attachments such as buccal tubes, sheaths, cleats and lugs.

**Description of the device**

Molar and Bicuspid bands are manufactured from annealed stainless steel shaped to fit the anatomy of the tooth, both maxillary and mandibular. They are permanently etch marked on the exterior with the size and quadrant. The interior is micro-etched for better retention.

**1. Part Numbers**

Product	Part	Variants	Material
Molar Bands	B10* K10* B11* K11* B12* K12* B13* K13		

**Warnings and Precautions**

	Devices are designed and manufactured for single patient use only. The reuse or reprocessing may cause cross infection, loss of material or mechanical integrity and strength.
	Avoid cementing bands in occlusion to prevent wearing of enamel surfaces during all stages of treatment due to the hardness of the bands.
	Cementing to compromised teeth can increase the risk of tooth damage.
	Damage to the enamel of a tooth or to a restoration (crown, bonding, veneer, etc.) is possible when orthodontic appliances are removed.
	Appliance breakage and or distortion requires immediate attention by the orthodontist.
	Patients should wear a mouthguard when engaging in sports activities; injury to the patient and /or breakage or distortion of the appliance can occur. In the event of injury or ingestion, instruct the patient to seek medical care.
	Care should be taken to avoid habits or foods that could break the bands. Follow orthodontist recommendations for proper at home care. In the event of injury or ingestion, instruct the patient to seek medical care.
	This product contains nickel or chromium, a chemical known to the State of California to cause cancer. In the event of an allergic reaction, immediately remove orthodontic appliance and instruct the patient to seek medical care.
	This medical device is MR Conditional. Scanning of the head and neck may cause image artifacts. If an MRI is needed, contact the radiologist prior to the procedure for instructions regarding appliance removal, if required.

**Clinical Benefits**

Bands, as a part of the orthodontic system, brings an aesthetic improvement with the alignment and leveling of the teeth, better hygiene, more efficient masticatory function facilitating swallowing, better digestion, prevention of occlusal problems, better speech and breathing.

**Residual Risks**

Orthodontic treatments are subject to relapse. Bands, as part of the fixed orthodontic system, can contribute to discoloration or decalcification of teeth, bone and root resorption, and periodontal complications

**Potential Side Effects or Adverse Events**

1. Allergic reaction, particularly if band contains nickel.
2. Discomfort and pain
3. Oral and mucosal damage
4. Ingestion
5. Aspiration
6. Adverse tooth movement resulting from wire distortion
7. Bands, as part of the fixed orthodontic system, can contribute to difficulties in maintaining oral hygiene, difficulty speaking or chewing and halitosis.

**Instructions for use**

The use of Orthodontic bands is restricted to dental professionals with qualification in orthodontics. Use by persons outside this context may result in irreversible damage to the patient. The choice of size and shape to be used is made by the qualified professional, based on their knowledge acquired during their training in conjunction with the diagnosis and treatment plan for each patient.

	To avoid ingestion or aspiration use of an appropriate barrier is recommended when seating, cementing, and removing bands.
---	--

1. Select the size of the band, choosing the one that best fits the tooth.
2. If welding a tube or attachment, weld onto the band according to knowledge acquired during training and coordinate with bracket placement.
3. Prepare the tooth surface according to the instructions of the manufacturer of the acquired orthodontic adhesive or cement.
4. Apply the resin or orthodontic cement inside the band, according to the manufacturer's instructions for use of the adhesive or cement purchased,
5. Place the band on the teeth using specific instruments, according to the manufacturer's instructions for use of the adhesive or cement purchased and knowledge acquired during training.

	When using a band pusher, exercise caution to prevent the instrument from slipping from the band thus injuring the patient's mucous membrane.
---	---

6. Remove excess cement from around the band and from the occlusal surface of the tooth.
7. After treatment, remove the band and orthodontic adhesive from the teeth using an appropriate instrument as per the orthodontic adhesive manufacturer's instructions for use and knowledge acquired during training.

Patient should be monitored regularly when undergoing orthodontic treatment. Particular care should be taken to monitor for band loosening or damage.

Any serious incident experienced with the use of the product should be reported to the manufacturer and the competent authority.

## Maintenance and Storage instructions

1. Store bands in product packaging until use.
2. The professional should use good clinical practices for the disposal of medical products in accordance with current sanitary regulations.

## Patient Instructions

	Patients should wear a mouthguard when engaging in sports activities; injury to the patient and /or breakage or distortion of the appliance can occur. In the event of injury or ingestion, instruct the patient to seek medical care.
	Follow orthodontist recommendations for proper at home care.
	Appliance breakage and or distortion requires immediate attention by the orthodontist.
	Care should be taken to avoid habits or foods that could break the appliance. In the event of injury or ingestion, instruct the patient to seek medical care.
	This product contains nickel or chromium, a chemical known to the State of California to cause cancer. In the event of an allergic reaction, immediately remove orthodontic appliance and instruct the patient to seek medical care.
	This medical device is MR Conditional. Scanning of the head and neck may cause image artifacts. If an MRI is needed, contact the radiologist prior to the procedure for instructions regarding appliance removal, if required.
	Practice good oral hygiene to prevent plaque installation and subsequent periodontal disease.

## Warranty

Limited Warranty. Product(s) manufactured by FlexMedics are warranted to be free of defects in materials and manufacture for ninety (90) days from date of delivery. Any material or manufacture defect covered by this limited warranty which occurs during normal use and is reported to Seller in writing during the period of ninety (90) days from the date such products are shipped to customer. FlexMedic's obligation hereunder, upon verification of the defect or error, shall be to provide one of the following, in FlexMedic's sole discretion: (i) replacement at no charge to customer; or (ii) credit the purchase price to customer. If, upon the inspection of FlexMedics product to which this warranty applies, FlexMedics determines that a claimed defect was not due to its manufacture or materials, FlexMedics will invoice customer for the value of product and services at customer's expense. This warranty shall not apply to any product that has been subject to misuse, misapplication, neglect (including but not limited to use of unauthorized parts or attachments), adjustment or repair and/or improper handling and storage. Identifiable items manufactured by others but installed in or affixed to FlexMedics products, or distributed by FlexMedics to customer, are not warranted by FlexMedics but bear only the express warranty, if any, of the manufacturers thereof. No representative of FlexMedics has any authority to waive, alter, vary or add to the terms hereof without prior approval in writing, to customer, signed by an officer of FlexMedics. FlexMedic's liability for its products, whether for breach of contract, negligence, strict liability in tort, or otherwise, shall be limited, in FlexMedic's sole discretion, to the replacement of the products or parts thereof or a credit to the purchase price to customer. FlexMedics will not be liable for any other injury, loss, damage or expense, whether direct or consequential, including but not limited to loss of use, income, profit or production, or increased cost of operation, or spoilage of or damage to material, arising in connection with the sale, installation, use of, inability to use, or the repair or replacement of, or late delivery of, FlexMedics products.

Limited Warranty. Product(s) manufactured by Highland Metals are warranted to be free from defects. Products must be returned within 90 days from the date of purchase and must be returned in original packaging and in resalable condition. HMI™ reserves the right to final determination as to acceptance of product exchange or credit. Restock fee may apply. Prior to the return of any merchandise, please

contact HMI™ Customer Service for a return authorization number.

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