


TitanMoly™ – Nickel-Free Titanium Molybdenum Wire


Precision engineered and manufactured by G&H® in the U.S.A.



TitanMoly™ Titanium Molybdenum



ANSI/ADA Specifications #32, Orthodontic Wires: 2017
Force values expressed at 3mm deflection.



G&H QUALITY MEANS:

- Smoother surface finish for improved sliding mechanics
- Precise tolerances for optimal torque expression
- Accepts adaptation bends without breakage
- Nickel-free for sensitive patients
- Provides nearly twice the working range of stainless steel with 45% less deflection force and superior recovery




Call **800-526-1026** to order
or email us at orders@ghortho.com

See and order our full line of products at GHOrthodontics.com

TitanMoly™ titanium molybdenum wire exhibits a superior working range for various treatment adaptations. All orders of treatment bends can be applied without breakage.

TitanMoly™ Titanium Molybdenum vs. Stainless Steel 3-Point Deflection Test



■ TitanMoly™ Titanium Molybdenum ■ Stainless Steel

.016X.022 wire sizes used in this comparison



GHOrthodontics.com

TOLL-FREE 800-526-1026

orders@ghortho.com

G&H Orthodontics 40 Linville Way Franklin, IN 46131 U.S.A. 800-526-1026 U.S. & Canada +1 317-346-6655 International

TitanMoly™ – Nickel-Free

Titanium Molybdenum Wire

Precision engineered and manufactured by G&H® in the U.S.A.



ORTHODONTICS®
Quality Manufacturing. Exceptional Service.



GHOrthodontics.com

TOLL-FREE 800-526-1026

orders@ghortho.com

G&H Orthodontics 40 Linville Way Franklin, IN 46131 U.S.A. 800-526-1026 U.S. & Canada +1 317-346-6655 International