

ARCHWIRES

G4™

NICKEL TITANIUM

Precision engineered and manufactured
by G&H® in the U.S.A.

G&H QUALITY MEANS:


- Low load forces of engagement
- Easy wire deflection and delivery to bracket
- Less stress on bonds
- Increased patient comfort
- High residual forces of recovery




TOLL-FREE 800-526-1026

orders@ghortho.com


G4™ NICKEL TITANIUM



Three-point deflection testing to measure force values



ANSI/ADA Specifications #32, Orthodontic Wires: 2017



Description: Nickel Titanium .016 – Test Temp: 36 °C – Super Elastic (SE/G4) fully activated by 25°C



ORTHODONTICS®
Quality Manufacturing. Exceptional Service.

TOLL-FREE 800-526-1026
orders@ghortho.com

GHorthodontics.com

G&H Orthodontics 40 Linville Way Franklin, IN 46131 U.S.A.
800-526-1026 U.S. & Canada +1 317-346-6655 International