

ARCHWIRES

M5™


THERMAL COPPER

NICKEL TITANIUM

Precision engineered and manufactured
by G&H® in the U.S.A.

G&H QUALITY MEANS:


- Gentle loading forces for easier engagement
- Precise temperature transformation of 28°C
- Consistent and predictable tooth moving forces
- Proprietary mechanical high polish technique for reduced surface friction
- US manufacturing quality



TOLL-FREE 800-526-1026

orders@ghortho.com

M5™ THERMAL COPPER NiTi




ANSI/ADA Specifications #32, Orthodontic Wires: 2017 - Force values expressed at 3mm deflection.

Recovery Force (gram-force)

M5™ Thermal Copper NiTi x Leading Competitor's Copper NiTi

Force measured from 3-pt bend test



Tested to ANSI/ADA Specification #32, Orthodontic Wires: 2006



TOLL-FREE 800-526-1026
orders@ghortho.com

GHorthodontics.com

G&H Orthodontics 40 Linville Way Franklin, IN 46131 U.S.A.
 800-526-1026 U.S. & Canada +1 317-346-6655 International